



# Cochlear™ Nucleus® CP802 Sound Processor

Part of the market-leading Nucleus® family  
of products from the global leader in  
implantable hearing solutions

*Hear now. And always*



Cochlear®

The Cochlear™ Nucleus® CP802 Sound Processor is designed to withstand the rough and tumble of everyday life – so it's tough enough to keep up with your child as they learn, play and grow.

### Durable design

- ✓ Made from **titanium and medical grade plastic**
- ✓ **Water resistant (IP57)**<sup>1</sup> and splashproof for peace of mind
- ✓ **Rated IP68** when worn with the Cochlear Aqua Accessory so children can enjoy their favourite water activities

### Everyday hearing

- ✓ **Two omni-directional microphones** to pick up sound no matter where it's coming from.
- ✓ Cochlear's SmartSound® Everyday program **optimises the sound signal** to help your child hear clearly and feel connected in everyday situations.
- ✓ Inbuilt telecoil **improves your child's access to sound** by allowing wireless connection to phones and places with loop systems (such as schools or cinemas).
- ✓ Audio input allows **easy connection to headphones, portable music systems or FM systems** used in schools – to help your child hear their teacher's voice more clearly.

### Parental control at the touch of a button

- ✓ Control your child's processor settings from the device itself or remotely using the **Cochlear Nucleus CR230 Remote Assistant** – so there's no need to interrupt what your child is doing.
- ✓ Multi-colour LED display on the Nucleus CP802 Sound Processor uses different lights to **signal to parents** that the processor is on and the microphone is receiving sound, or that the processor needs attention – for example, if it has a flat battery.
- ✓ **Tamper-proof processor buttons** are lockable to prevent children from accidentally changing settings.



The durable  
Cochlear Nucleus CP802  
Sound Processor




The convenient  
Cochlear Nucleus CR230  
Remote Assistant

### Your trusted partner

When you choose Cochlear's Nucleus CP802 Sound Processor, you are choosing the global leader in implantable hearing solutions.

We've helped over 450,000 people just like you, who have chosen Cochlear's solutions. So you can feel confident that wherever you are in the world or on your hearing journey, Cochlear is here to support you.



# Hear now. And always

As the global leader in implantable hearing solutions, Cochlear is dedicated to bringing the gift of sound to people with moderate to profound hearing loss. We have helped over 450,000 people of all ages live full and active lives by reconnecting them with family, friends and community.

We aim to give our recipients the best lifelong hearing experience and access to innovative future technologies. For our professional partners, we offer the industry's largest clinical, research and support networks.

That's why more people choose Cochlear than any other hearing implant company.

## Reference

1. Cochlear Limited. Nucleus CP802 Sound Processor User Guide. Data on file.

 **Cochlear Ltd** (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia T: +61 2 9428 6555 F: +61 2 9428 6352

[www.cochlear.com](http://www.cochlear.com)

Please seek advice from your medical practitioner or health professional about treatments for hearing loss. They will be able to advise on a suitable solution for the hearing loss condition. All products should be used only as directed by your medical practitioner or health professional. Outcomes and results may vary. Not all products are available in all countries. Please contact your local Cochlear representative.

ACE, Advance Off-Stylet, AOS, AutoNRT, Autosensitivity, Beam, Button, CareYourWay, Carina, Cochlear, コクレア, 科利耳, Cochlear SoftWear, Codacs, ConnectYourWay, Contour, Contour Advance, Custom Sound, ESPrit, Freedom, Hear now. And always, HearYourWay, Hugfit, Hybrid, Invisible Hearing, Kanso, MET, MicroDrive, MP3000, myCochlear, mySmartSound, NRT, Nucleus, Off-Stylet, Slimline, SmartSound, Softip, SPrint, True Wireless, the elliptical logo, WearYourWay and Whisper are either trademarks or registered trademarks of Cochlear Limited. Ardium, Baha, Baha SoftWear, BCDrive, DermaLock, EveryWear, Vistafix, and WindShield are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB.

© Cochlear Limited 2018. D1517038 ISS1 OCT18

