

Live life your way

Choose your Cochlear™ sound processor

Hear now. And always



Only Cochlear sound processors let you **RECAPTURE** the joy of **A+ HEARING**

Kanso™ is Cochlear's newest sound processor. It's discreet and easy to use, without compromising on hearing performance technology.

Cochlear offers a choice of two sound processors in the Nucleus portfolio :

- Kanso
- Nucleus 6

These processors share many of the same features and benefits, but they also have some important differences, giving each processor distinct advantages. Kanso is designed to maximise discretion and comfort, while Nucleus 6 is designed for a feature-rich user experience.



KANSO™

A smart, simple, discreet, off-the-ear sound processor with all the proven technology of the Nucleus 6 Sound Processor.



NUCLEUS® 6

The Nucleus 6 System is our most advanced system ever, featuring breakthrough technologies designed to make your hearing journey to A+ hearing a smooth experience.

Leap into the world of **HEARING WITH EASE**

SmartSound® iQ processing technology replicates the functions of a normal hearing brain, so you always automatically receive crisp and clear sounds, wherever you are. First the dual microphones detect the sounds and then the processor cleans and adjusts them for optimal clarity, without you having to do a thing. That means you can go from the noisy street to the quiet office without having to think about your hearing. Just like everyone else.



Leap into the world of **SUPER HEARING**

Cochlear™ Wireless Accessories stream directly to your Nucleus 6 and Kanso Sound Processor so you can hear clearly in situations where even people with normal hearing might struggle. You can take a phone call in a busy train station or stream a lecturer's voice straight to your sound processor – cutting through the background noise.




 KANSO™	BOTH FEATURE:	 NUCLEUS® 6
Off-the-ear sound processor	 SMARTSOUND IQ WITH SCAN	Behind-the-ear sound processor
Lightest and smallest off-the-ear sound processor	 DUAL MICROPHONES	Lightest and smallest behind-the-ear sound processor
Eight colours 	 DUST AND SPLASH PROOF	Five colours 
Disposable batteries	 WATERPROOF WITH AQUA+	Disposable and rechargeable batteries
No acoustic component	 COMPATIBLE WITH TRUE WIRELESS™ DEVICES	Offers Hybrid™ Hearing with acoustic component
Not compatible with Nucleus 22 implants	 CR210 REMOTE CONTROL AND CR230 REMOTE ASSISTANT	Compatible with all Nucleus implants

Whichever device you choose, you can be assured that you're getting the great hearing performance of SmartSound® iQ with SCAN, dual microphone directionality and True Wireless™ technology.

The Cochlear difference: A+ HEARING

A+ HEARING is only possible when you combine the four pillars of hearing performance; perfected over our many years as the global market leader. This knowledge is further strengthened by the experience shared by more than 450,000 people we have treated since inventing the world's first cochlear implant nearly 40 years ago.

FOUR PILLARS OF HEARING PERFORMANCE



Hear now. And always

As the global leader in implantable hearing solutions, Cochlear is dedicated to bringing the gift of sound to people with moderate to profound hearing loss. We have helped over 450,000 people of all ages live full and active lives by reconnecting them with family, friends and community.

We give our recipients the best lifelong hearing experience and access to innovative future technologies. For our professional partners, we offer the industry's largest clinical, research and support networks.

That's why more people choose Cochlear than any other hearing implant company.

Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia Tel: +61 2 9428 6555 Fax: +61 2 9428 6352.
Cochlear NZ Limited Level 4, Takapuna Towers, 19-21 Como Street, Takapuna, Auckland 0622 New Zealand Tel: 64 9 914 1983 Fax: 0800 886 036.
Nihon Cochlear Co Ltd Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 0245.
Cochlear (HK) Ltd Room 1204, 12/F, CRE Building, No 303 Hennessey Road, Wanchai, Hong Kong SAR Tel: +852 2530 5773 Fax: +852 2530 5183.
Cochlear Medical Device (Beijing) Co Ltd Unit 2208-2212, Tower B, Gemdale Building No.91 Jianguo Road, Chaoyang District Beijing 100022 P.R. CHINA Tel: +86 10 5909 7800 Fax: +86 10 5909 7900.
Cochlear Limited (Singapore Branch) 238A Thomson Road, #25-06, Novena Square Office Tower A, Singapore 307684 Tel: +65 6553 3814.
Cochlear Korea Ltd 1st floor, Cheongwon building, 828-5, Yuksam dong, Kangnam gu, Seoul, Korea Tel: +82 2 533 4450 Fax: +82 2 533 8408.
Cochlear Medical Device Company India Private Ltd Ground Floor, Platina Building, Plot No C-59, G-Block, Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India Tel: +91 22 6112 1111 Fax: +91 22 6112 1100.

www.cochlear.com

Connect with us on Social media:  /CochlearAsiaPacific  @CochlearAsiaPac  @CochlearAsiaPacific

Nucleus and Hear now. And always, Cochlear and the elliptical logo are trademarks or registered trademarks of Cochlear Limited. © Cochlear Limited 2016. D846281 ISS1 JUN16

