

KANSO™

Product Catalogue



KANSO™



Cochlear®

Hear now. And always

Introducing **KANSO**, the discreet
new sound processor from Cochlear

KANSO

Kanso™ is an off-the-ear sound processor that's so discreet it's barely noticeable, so small and comfortable you can forget you're wearing it.

This offers greater confidence that you can live your life your way.

Kanso does this without compromising on hearing performance technology. Kanso incorporates the same industry-leading features as the Cochlear™ Nucleus® 6 Sound Processor, making it even easier for you to find the joy of hearing.



DISCREET

It's small and worn off-the-ear

SIMPLE

Easy to use, you can set and forget

SMART

SmartSound® iQ with SCAN
and dual microphones

To order items in this product catalogue
please contact your clinic or nearest
Cochlear office.

Please note that:

- Images shown are not actual size.
- Products are available as single items only, unless otherwise stated.
- Products shown in this catalogue may not be available in all countries and configurations or product specifications may vary.

Only your expression gives Kanso away

KANSO (CP950) SOUND PROCESSOR

Kanso is available in eight colours to blend in with your hair.

Part Number:

P771862
White

P771872
Silver

P771867
Slate Grey

P771866
Black



P771870
Sandy Blonde

P771863
Golden Blonde

P771875
Copper Brown

P771873
Chocolate Brown



MAGNETS

There are seven magnet strengths to choose from to suit your needs.

Part Number:

Z566412 CP950 Magnet 1/2M

Z502922 CP950 Magnet 1M

Z502923 CP950 Magnet 2M

Z502924 CP950 Magnet 3M

Z502925 CP950 Magnet 4M

Z566414 CP950 Magnet 5M

Z566415 CP950 Magnet 6M



REMOTE CONTROL AND REMOTE ASSISTANT

Adjust simple settings without having to touch your Kanso thanks to these discreet, easy to use wireless devices.

CR210 Remote Control

A pocket-sized remote control with a handy display screen for changing simple settings.

Part Number:

Z418890 CR210 Remote Control



CR230 Remote Assistant

The intuitive design of this full-featured device puts you in control of monitoring and managing your own hearing.

Part Number:

Z418891 CR230 Remote Assistant



REMOTE ACCESSORIES

CR200 Series Remote Assistant Charging Kit

Part Number:

Z319245 CR200 Series Remote Assistant Charging Kit



CR200 Series Protective Case (CR230)

Supplied as a pack of 2, includes a lanyard.

Part Number:

Z319244 CR200 Series Protective Case (CR230)



Remote Assistant Belt Case

Includes lanyard.

Part Number:

Z332268 Remote Assistant Belt Case



CR200 Series USB Cable (Micro)

Replacement cable for CR200 Series Charging Kit.

Part Number:

Z319241 CR200 Series USB Cable (Micro)



ACTIVE LIFESTYLES

Available in different sizes and with hair or clothing clips for added retention.

Nucleus Safety Line - Standard

Part Number:		Length
Z467062	Nucleus Safety Line - Standard (2 pack)	21 cm



Nucleus Safety Line - Long

Part Number:		Length
P742062	Nucleus Safety Line - Long (2 pack)	26 cm



Nucleus Safety Line - Short Double Loop


Part Number:		Length
P798366	Nucleus Safety Line - Short Double Loop (2 pack)	8 cm



Nucleus Safety Line - Short Double Loop with Hair Clips:

Includes 2 Safety Lines and 5 Hair Clips.

Part Number:		Length	Hair Clip Colour
P743011	Nucleus Safety Line - Short Double Loop	8 cm	Black
P789715	Nucleus Safety Line - Short Double Loop	8 cm	Brown
P789713	Nucleus Safety Line - Short Double Loop	8 cm	White



COCHLEAR HEADBAND

Choose one of six headbands to suit your lifestyle.

Part Number:		Size	Colour	Average head circumference
P705126	Cochlear Headband	XXS	Vanilla	41 cm
P783375	Cochlear Headband	XS	Pink	47 cm
P783380	Cochlear Headband	S	Dark Blue	50 cm
P783385	Cochlear Headband	M	Purple	54 cm
P783387	Cochlear Headband	M	Black	54 cm
P783388	Cochlear Headband	L	Black	58 cm



CARE AND COMFORT

SoftWear Pads

The SoftWear Pads cushion the Kanso Sound Processor against the head for extra comfort.

Part Number:
P793406 SoftWear Pads (20 pads)



Microphone Protector Kit

We recommend changing the microphone protectors every three months to keep hearing at your best.

Part Number:
P740443 Microphone Protector Kit (2 pack)
P743057 Microphone Protector Kit (4 pack)



Storage Case

The Storage Case can hold two processors.

Part Number:
Z342010 Storage Case



Battery Case Keyring

This carry case holds your spare batteries.

Part Number:
Z322700 Battery Case Keyring holds 2 batteries



Breeze by Dry & Store®

This user-friendly product removes moisture from your sound processor. It helps to increase the processor's longevity and may result in fewer repairs and better sound quality. Available in regional versions to suit specific power supply and voltage requirements.*

Part Number:
Z60761 Breeze by Dry & Store (Australia/New Zealand)
Z267427 Breeze by Dry & Store (Korea)
Z295319 Breeze by Dry & Store (Japan)



Dry-Brik® by Dry & Store

This product is to be used in conjunction with the Dry & Store unit to remove moisture from your sound processor.

Part Number:
Z60764 Dry-Brik (3 pack)



Drying Capsules

For use in your everyday Storage Case to remove moisture from your sound processor, the drying capsule changes colour when it needs to be replaced.

Part Number:
N94527 Drying Capsules (4 pack)



Metal Can Desiccant

The Metal Can Desiccant works similarly to the Dry-Brik Desiccant, absorbing moisture during the drying process. The Metal Can Desiccant is reusable.

Part Number:
Z61088 Metal Can Desiccant

BATTERY INFORMATION

Kanso uses disposable batteries that last for up to 60 hours.* Batteries are supplied as a pack of six batteries on one card.

Cards can be sold separately or in boxes of 10 cards (60 batteries) or in a carton of 5 boxes (300 batteries).

Part Number:
B13033 Power One Implant Plus P675



*Supplied with global power adaptor and one Dry-Brik® Desiccant (Z60764).

*Battery life will vary depending on the power requirements of your listening program and your individual characteristics (e.g. use of accessories.)

ACTIVITY KIT

AQUA+ FOR KANSO

The Aqua+ is a soft, flexible silicone sleeve that fits over Kanso Sound Processor, making it completely waterproof.* It has an ingress protection rating of IP68.

CP950 Aqua+ for Kanso

Part Number:
Z594543 CP950 Aqua+ for Kanso (2 pack)



Battery LR44

Part Number:
P741235 Battery LR44 (pack of 2 batteries)



Activity Kit Case

Part Number:
P740569 Activity Kit case



Nucleus Safety Line - Long

Part Number:	Length
P742062 Nucleus Safety Line - Long (2 pack)	26 cm



Cochlear Headband

Choose one of six headbands to suit your lifestyle.

Part Number:	Size	Colour	Average head circumference
P705126 Cochlear Headband	XXS	Vanilla	41 cm
P783375 Cochlear Headband	XS	Pink	47 cm
P783380 Cochlear Headband	S	Dark Blue	50 cm
P783385 Cochlear Headband	M	Purple	54 cm
P783387 Cochlear Headband	M	Black	54 cm
P783388 Cochlear Headband	L	Black	58 cm



COCHLEAR™ TRUE WIRELESS™ RANGE

Stream conversation, music and TV direct to your Kanso Sound Processor, providing the best possible hearing even in noise and over distance.

Phone Clip

The Phone Clip gives you a hands-free connection to friends and family using your smartphone. Make calls, carry out conversations, stream music and adjust volume at the touch of a button.

The product will be delivered in a kit that includes the accessory with: USB cable (USB – micro USB), AC power supply, pouch and lanyard.

Part Number:

94770	Phone Clip (Europe)
94771	Phone Clip (US)
94772	Phone Clip (Great Britain)
94773	Phone Clip (Australia/ New Zealand)



Mini Microphone 2+

The Mini Microphone 2+ transmits speech and audio directly to sound processors across a range of up to 25 metres (clear line of sight). So whether you're relaxing at home or getting out and being active, you can be closer to the sounds, experiences and people that matter most.

The product will be delivered in a kit that includes the accessory with: Charger adaptor, charging cord (micro-USB), jack cable, pouch, lanyard.

Part Number:

P770846	Mini Microphone 2+ (Europe)
P770848	Mini Microphone 2+ (US)
P770845	Mini Microphone 2+ (Great Britain)
P770847	Mini Microphone 2+ (Australia/ New Zealand)



TV Streamer

Enjoy crisp, clear stereo sound directly from your TV. With the TV Streamer, you can enjoy your favourite show while still being part of any conversation going on around you.

The product will be delivered in a kit that includes the accessory with: AC/DC adaptor, USB jack, scart, phone adaptor, USB cable (USB-micro USB) and jack-phone cable red-white.

Part Number:

94760	TV Streamer (Europe)
94761	TV Streamer (US)
94762	TV Streamer (Great Britain)
94763	TV Streamer (Australia/ New Zealand)



*You must use alkaline, nickel metal hydride or silver oxide disposable batteries with your Aqua+ for Kanso.

Hear now. And always

As the global leader in implantable hearing solutions, Cochlear is dedicated to bringing the gift of sound to people with moderate to profound hearing loss. We have helped over 450,000 people of all ages live full and active lives by reconnecting them with family, friends and community.

We aim to give our recipients the best lifelong hearing experience and access to innovative future technologies. For our professional partners, we offer the industry's largest clinical, research and support networks.

That's why more people choose Cochlear than any other hearing implant company.

 Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

www.cochlear.com

Follow us on:



Please seek advice from your medical practitioner or health professional about treatments for hearing loss. They will be able to advise on a suitable solution for the hearing loss condition. All products should be used only as directed by your medical practitioner or health professional. Not all products are available in all countries. Please contact your local Cochlear representative.

ACE, Advance Off-Stylet, AOS, AutoNRT, Autosensitivity, Beam, Button, Carina, Cochlear, コクレア, Codacs, Contour, Contour Advance, Custom Sound, ESPrit, Freedom, Hear now. And always, Hugfit, Hybrid, inHear, Invisible Hearing, Kanso, MET, MP3000, myCochlear, NRT, Nucleus, 科利耳, Off-Stylet, SmartSound, Softip, SPrint, the elliptical logo, True Wireless and Whisper are either trademarks or registered trademarks of Cochlear Limited. Ardium, Baha, Baha Divino, Baha Intenso, Baha PureSound, Baha SoftWear, BCDrive, DermaLock, Vistafix and WindShield are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB.

© Cochlear Limited 2016.
DB84576 ISS3 JULY16

